


HANAROSTENT® Biliary Lasso (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the both ends for Anti-Migration
- Fully covered stent
- Radiopaque markers for fluoroscopic guidance
- Lasso on the proximal for stent removal and long lasso for Inside-Out stent removal

Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCBC08040-B180 ⋮ BCBC08120-B180	11.5-8-11.5	26~106*	40~120*	1800	10.2
BCBC10040-B180 ⋮ BCBC10120-B180	13.5-10-13.5	26~106*	40~120*	1800	10.2

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device

Distal Proximal


(NNN) Non covered
(CCC) Fully covered
(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length
② Stent Diameter ④ Delivery Device Length

HANAROSTENT® Biliary Lasso (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the both ends for Anti-Migration
- Fully covered stent
- Radiopaque markers for fluoroscopic guidance
- Short lasso on the proximal for stent removal and long lasso on the distal for Inside-Out stent removal

Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCBC08040-N060 ⋮ BCBC08120-N060	11.5-8-11.5	26~106*	40~120*	600	10.2
BCBC10040-N060 ⋮ BCBC10120-N060	13.5-10-13.5	26~106*	40~120*	600	10.2

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length