


HANAROSTENT® Biliary Flap (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- 4 Anchoring Flaps on the both ends for Anti-Migration
- Radiopaque markers for fluoroscopic guidance
- Fully covered stent


Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCG-06-040-180 ⋮ BCG-06-120-180	6	20~100*	40~120*	1800	8
BCG-08-040-180 ⋮ BCG-08-120-180	8	20~100*	40~120*	1800	8.5
BCG-10-040-180 ⋮ BCG-10-120-180	10	20~100*	40~120*	1800	8.5

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④


① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length

HANAROSTENT® Biliary Flap (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- 4 Anchoring Flaps on the both ends for Anti-Migration
- Radiopaque markers for fluoroscopic guidance
- Fully covered stent

Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCG-06-040-060 ⋮ BCG-06-120-060	6	20~100*	40~120*	600	8
BCG-08-040-060 ⋮ BCG-08-120-060	8	20~100*	40~120*	600	8.5
BCG-10-040-060 ⋮ BCG-10-120-060	10	20~100*	40~120*	600	8.5

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length