


HANAROSTENT® Biliary Flap (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the proximal end for Anti-Migration
- 4 Anchoring Flaps on the distal part for Anti-Migration
- Radiopaque markers for fluoroscopic guidance
- Fully covered stent

Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCS08040-D180 ⋮ BCS08120-D180	11.5-8	25-105*	40~120*	1800	8.5
BCS10040-D180 ⋮ BCS10120-D180	13.5-10	25-105*	40~120*	1800	8.5

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device

Distal Proximal


(NNN) Non covered
(CCC) Fully covered
(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length
② Stent Diameter ④ Delivery Device Length

HANAROSTENT® Biliary Flap (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the distal end for Anti-Migration
- 4 Anchoring Flaps on the proximal part for Anti-Migration
- Radiopaque markers for fluoroscopic guidance
- Fully covered stent

Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCS08040-R060 ⋮ BCS08120-R060	8-11.5	25-105*	40~120*	600	8.5
BCS10040-R060 ⋮ BCS10120-R060	10-13.5	25-105*	40~120*	600	8.5

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length