


HANAROSTENT® Biliary Flap Lasso (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the proximal and 4 Anchoring Flaps on distal end for Anti-Migration
- Fully covered stent
- Radiopaque markers for fluoroscopic guidance
- Lasso on the proximal for stent removal and long lasso for Inside-Out stent removal

Endoscopic Application


	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCT-08-040-180 ⋮ BCT-08-120-180	8-11.5	23~103*	40~120*	1800	9
BCT-10-040-180 ⋮ BCT-10-120-180	10-13.5	23~103*	40~120*	1800	9

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

<p>• Based on Delivery Device</p>		<p>(NNN) Non covered</p> <p>(CCC) Fully covered</p> <p>(NCN) Partially covered</p>	<p>XXX-00-000-000 XXX 00 000-X000</p> <p>① ② ③ ④ ① ② ③ ④</p> <p>① Product Name ③ Stent Length</p> <p>② Stent Diameter ④ Delivery Device Length</p>
-----------------------------------	--	--	--

HANAROSTENT® Biliary Flap Lasso (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the distal and 4 Anchoring Flaps on the proximal end for Anti-Migration
- Fully covered stent
- Radiopaque markers for fluoroscopic guidance
- Short lasso on the proximal for stent removal and long lasso on the distal for Inside-Out stent removal

Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCT-08-040-060 ⋮ BCT-08-120-060	8.11.5	23~103*	40~120*	600	9
BCT-10-040-060 ⋮ BCT-10-120-060	10-13.5	23~103*	40~120*	600	9

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length