


HANAROSTENT® Multi-Hole™ Biliary (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the both ends for Anti-Migration
- Multi holes on the membrane to facilitate bile juice drainage
- Radiopaque markers for fluoroscopic guidance
- Lasso on the proximal for stent removal

Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BMCL-08-040-180 ⋮ BMCL-08-120-180	11.5-8-11.5	26~106*	40~120*	1800	8
BMCL-10-040-180 ⋮ BMCL-10-120-180	13.5-10-13.5	26~106*	40~120*	1800	8

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④


① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length

HANAROSTENT® Multi-Hole™ Biliary (CCC)



Indication :

- Palliative treatment of bile duct malignant and/or benign strictures

Features :

- Flare on the both ends for Anti-Migration
- Multi holes on the membrane to facilitate bile juice drainage
- Radiopaque markers for fluoroscopic guidance
- Lasso on the proximal for stent removal


Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BMCL-08-040-060 ⋮ BMCL-08-120-060	11.5-8-11.5	26~106*	40~120*	600	8
BMCL-10-040-060 ⋮ BMCL-10-120-060	13.5-10-13.5	26~106*	40~120*	600	8

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length