


HANAROSTENT® Biliary Valve (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Flare on the distal end for Anti-Migration
- Uncovered stent
- Radiopaque markers for fluoroscopic guidance
- S-Valve for Anti-Reflux

Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BNX-08-060-180 ⋮ BNX-08-120-180	11.5-8	53~113*	60~120*	1800	10.2
BNX-10-060-180 ⋮ BNX-10-120-180	13.5-10	53~113*	60~120*	1800	10.2

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~120 : 60,70,80,90,100,110,120)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④


① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length

HANAROSTENT® Biliary Valve (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Flare on the proximal end for Anti-Migration
- Uncovered stent
- Radiopaque markers for fluoroscopic guidance
- S-Valve for Anti-Reflux


Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BNX-08-060-060 ⋮ BNX-08-120-060	8-11.5	53~113*	60~120*	600	10.2
BNX-10-060-060 ⋮ BNX-10-120-060	10-13.5	53~113*	60~120*	600	10.2

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~120 : 60,70,80,90,100,110,120)

• Based on Delivery Device



Distal

Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length