


## HANAROSTENT® Biliary (NCN)



### Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

### Features :

- Flare on the both ends for Anti-Migration
- Uncovered on the both ends for Anti-Migration
- Radiopaque markers for fluoroscopic guidance

### Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BPA08040-D180 ⋮ BPA08120-D180	11.5-8-11.5	26~106*	40~120*	1800	8.5
BPA10040-D180 ⋮ BPA10120-D180	13.5-10-13.5	26~106*	40~120*	1800	8.5

※ Sizing and availability varies by country

\* Increment of 10mm ( Total length 40~120 : 40,50,60,70,80,90,100,110,120 )

• Based on Delivery Device

Distal Proximal


(NNN) Non covered  
(CCC) Fully covered  
(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length  
② Stent Diameter ④ Delivery Device Length

## HANAROSTENT® Biliary (NCN)



### Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

### Features :

- Flare on the both ends for Anti-Migration
- Uncovered on the both ends for Anti-Migration
- Radiopaque markers for fluoroscopic guidance


### Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BPA08040-R060 ⋮ BPA08120-R060	11.5-8-11.5	26~106*	40~120*	600	8.5
BPA10040-R060 ⋮ BPA10120-R060	13.5-10-13.5	26~106*	40~120*	600	8.5

※ Sizing and availability varies by country

\* Increment of 10mm ( Total length 40~120 : 40,50,60,70,80,90,100,110,120 )

• Based on Delivery Device



Distal

Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length