


CHOOSTENT® Esophagus ST (CCN)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :

- Dumbbell heads on the both ends for Anti-Migration
- Shim's Technique on the proximal end for Anti-Migration
- Uncovered on the proximal head for Anti-Migration
- To be positioned in the middle part of Esophagus
- Segmented body structure to conform to peristaltic movement
- Lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
EPG18060-X070 ⋮ EPG18170-X070	24-18-24	20~130*	60~170*	700	6/18
EPG22060-Z070 ⋮ EPG22170-Z070	28-22-28	20~130*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length