


CHOOSTENT® Esophagus Valve ST (CCN)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :


- Dumbbell heads on the both ends for Anti-Migration
- Uncovered on the proximal head for Anti-Migration
- S-Valve inside the stent for Anti-Reflux
- To be positioned in the Lower Esophageal Sphincter
- Shim's Technique on the proximal end for Anti-Migration
- Segmented body structure to conform to peristaltic movement
- Lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
HEVT-18-060-070 ⋮ HEVT-18-170-070	24-18-24	20~130*	60~170*	700	8/24
HEVT-22-060-070 ⋮ HEVT-22-170-070	28-22-28	20~130*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length