


HANAROSTENT® **Esophagus Valve** **BS (CCN)**



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :

- Bone shaped head on the proximal to minimize tissue overgrowth
- Uncovered on the proximal for Anti-Migration
- S-Valve inside the stent for Anti-Reflux
- To be positioned in the Lower Esophageal Sphincter
- Lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
EPV-18-060-070 ⋮ EPV-18-170-070	24-18-24	20~130*	60~170*	700	8/24
EPV-22-060-070 ⋮ EPV-22-170-070	28-22-28	20~130*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length