


CHOOSTENT® Esophagus (NCN)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :

- Dumbbell heads on the both ends for Anti-Migration
- Uncovered on the both ends for Anti-Migration
- To be positioned in the middle part of Esophagus
- Segmented body structure to conform to peristaltic movement
- Lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
EPE18060-X070 ⋮ EPE18170-X070	24-18-24	20~130*	60~170*	700	6/18
EPE20060-X070 ⋮ EPE20170-X070	26-20-26	20~130*	60~170*	700	6/18
EPE22060-Z070 ⋮ EPE22170-Z070	28-22-28	20~130*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length