


## CHOOSTENT® Esophagus Valve (CCN)



### Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

### Features :

- Dumbbell heads on the both ends for Anti-Migration
- S-Valve inside the stent for Anti-Reflux
- Uncovered on the proximal head for Anti-Migration
- To be positioned in the Lower Esophageal Sphincter
- Segmented body structure to conform to peristaltic movement
- Lassos on the both ends

|                                     | Stent(mm) |                |               | Delivery device |                 |
|-------------------------------------|-----------|----------------|---------------|-----------------|-----------------|
|                                     | Diameter  | Usable length* | Total length* | Length(mm)      | Diameter(mm/Fr) |
| EPF18060-Z070<br>⋮<br>EPF18170-Z070 | 24-18-24  | 20~130*        | 60~170*       | 700             | 8/24            |
| EPF22060-Z070<br>⋮<br>EPF22170-Z070 | 28-22-28  | 20~130*        | 60~170*       | 700             | 8/24            |

※ Sizing and availability varies by country

\* Increment of 10mm ( Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170 )

• Based on Delivery Device



Distal

Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length