


HANAROSTENT® Esophagus TTS (NCN)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :

- Dumbbell heads on the both ends for Anti-Migration
- To be positioned in the middle part of Esophagus
- Uncovered on the both ends for Anti-Migration
- Through The Scope (TTS)
- Lasso on the proximal end

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
EPBA-18-060-180 ⋮ EPBA-18-150-180	26-18-26	30~120*	60~150*	1800	3.5/10.5
EPBA-20-060-180 ⋮ EPBA-20-150-180	26-20-26	30~120*	60~150*	1800	3.5/10.5
EPBA-18-060-230 ⋮ EPBA-18-150-230	26-18-26	30~120*	60~150*	2300	3.5/10.5
EPBA-20-060-230 ⋮ EPBA-20-150-230	26-20-26	30~120*	60~150*	2300	3.5/10.5

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~150 : 60,70,80,90,100,110,120,130,140,150)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length