


HANAROSTENT® Esophagus Benign BS (CCC)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant and/or benign strictures

Features :


- Bone shaped head on the proximal to minimize tissue overgrowth
- To be positioned in the middle part of Esophagus
- Fully covered and lassos on the both ends
- Membrane outer covering for Easy Removal

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
EBN18060-Z070 ⋮ EBN18120-Z070	24-18-24	20~80*	60~120*	700	8/24
EBN20060-Z070 ⋮ EBN20120-Z070	26-20-26	20~80*	60~120*	700	8/24
EBN22060-Z070 ⋮ EBN22120-Z070	28-22-28	20~80*	60~120*	700	8/24

※ Sizing and availability varies by country

* Increment of 20mm (Total length 60~120 : 60,80,100,120)

• Based on Delivery Device



Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length