


HANAROSTENT®

Esophagus Benign

BS ST (CCC)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant and/or benign strictures

Features :

- Bone shaped head on the proximal to minimize tissue overgrowth
- Shim's Technique on the proximal end for Anti-Migration
- To be positioned in the middle part of Esophagus
- Fully covered and lassos on the both ends
- Membrane outer covering for Easy Removal

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
EBNT18060-Z070 ⋮ EBNT18120-Z070	24-18-24	20~80*	60~120*	700	8/24
EBNT20060-Z070 ⋮ EBNT20120-Z070	26-20-26	20~80*	60~120*	700	8/24
EBNT22060-Z070 ⋮ EBNT22120-Z070	28-22-28	20~80*	60~120*	700	8/24

※ Sizing and availability varies by country

*Increment of 20mm

* Increment of 20mm (Total length 60~120 : 60,80,100,120)

• Based on Delivery Device		(NNN) Non covered (CCC) Fully covered (NCN) Partially covered	XXX-00-000-000 ① ② ③ ④ ① Product Name ③ Stent Length ② Stent Diameter ④ Delivery Device Length	XXX 00 000-X000 ① ② ③ ④
-----------------------------------	--	--	--	---------------------------------------