


HANAROSTENT® **Esophagus Skidproof (CCC)**



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :

- Dumbbell heads on the both ends for anti-migration
- To be positioned in the middle part of Esophagus
- Fully covered and lassos on the both ends
- Skidproof body structure for Anti-Migration

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
ECD-20-080-070 ⋮ ECD-20-160-070	26-20-26	40~120*	80~160*	700	6/18

※ Sizing and availability varies by country

* Increment of 20mm (Total length 80~160 : 80,100,120,140,160)

<p>• Based on Delivery Device</p> <p>Distal Proximal</p>		<p>(NNN) Non covered</p> <p>(CCC) Fully covered</p> <p>(NCN) Partially covered</p>	<p>XXX-00-000-000</p> <p>① ② ③ ④</p> <p>① Product Name ③ Stent Length</p> <p>② Stent Diameter ④ Delivery Device Length</p>	<p>XXX 00 000-X000</p> <p>① ② ③ ④</p>
--	--	--	--	---------------------------------------