


## HANAROSTENT® Esophagus Lasso (NNN)



### Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

### Features :

- Uncovered stent
- Dumbbell heads on the both ends for Anti-Migration
- Lasso on the proximal end

|                                     | Stent(mm) |                |               | Delivery device |                 |
|-------------------------------------|-----------|----------------|---------------|-----------------|-----------------|
|                                     | Diameter  | Usable length* | Total length* | Length(mm)      | Diameter(mm/Fr) |
| END18060-C070<br>⋮<br>END18170-C070 | 24-18-24  | 20~130*        | 60~170*       | 700             | 4/12            |
| END20060-C070<br>⋮<br>END20170-C070 | 26-20-26  | 20~130*        | 60~170*       | 700             | 4/12            |

※ Sizing and availability varies by country

\* Increment of 10mm ( Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170 )

|                            |  |                         |                  |                          |
|----------------------------|--|-------------------------|------------------|--------------------------|
| • Based on Delivery Device |  | (NNN) Non covered       | XXX-00-000-000   | XXX 00 000-X000          |
|                            |  | (CCC) Fully covered     | ① ② ③ ④          | ① ② ③ ④                  |
|                            |  | (NCN) Partially covered | ① Product Name   | ③ Stent Length           |
|                            |  |                         | ② Stent Diameter | ④ Delivery Device Length |