


HANAROSTENT® Esophagus Bariatric Surgery ST (CCC)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant and/or benign strictures


Features :

- Shim's Technique on the proximal end for Anti-Migration
- Flexible zone in the middle part of stent for high conformity
- Fully covered and lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
ECBBT-28-180-090	32-28-28	145	180	900	8/24
ECBBT-28-210-090		175	210		
ECBBT-28-240-090		205	240		

※ Sizing and availability varies by country

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length

② Stent Diameter ④ Delivery Device Length