


HANAROSTENT® Esophagus Lasso (NNN)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant and/or benign strictures

Features :

- Large upper head for Anti-Migration
- Shim's Technique on the proximal end for Anti-Migration
- 2 Flexible zones in the middle part of stent for high conformity
- Round ending on the distal for less trauma
- Fully covered and lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
ECBST-26-180-090	36-26-34	130	180	900	8/24
ECBST-26-210-090		160	210		
ECBST-26-240-090		190	240		

※ Sizing and availability varies by country

• Based on Delivery Device



Distal

Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length