


HANAROSTENT® Biliary (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Flare on the both ends for Anti-Migration
- Uncovered stent
- Radiopaque markers for fluoroscopic guidance
- Large Cell body structure to facilitate second stent placement through the mesh


Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BNL-08-040-180 ⋮ BNL-08-100-180	11.5-8-11.5	26~86*	40~100*	1800	7
BNL-10-040-180 ⋮ BNL-10-100-180	13.5-10-13.5	26~86*	40~100*	1800	7

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~100 : 40,50,60,70,80,90,100)

• Based on Delivery Device



Distal

Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④


① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length

HANAROSTENT® Biliary (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Flare on the both ends for Anti-Migration
- Uncovered stent
- Radiopaque markers for fluoroscopic guidance
- Large Cell body structure to facilitate second stent placement through the mesh


Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BNL-08-040-060 ⋮ BNL-08-100-060	11.5-8-11.5	26~86*	40~100*	600	7
BNL-10-040-060 ⋮ BNL-10-100-060	13.5-10-13.5	26~86*	40~100*	600	7

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~100 : 40,50,60,70,80,90,100)

• Based on Delivery Device



Distal

Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length