

HANAROSTENT® Benefit™ Biliary (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Uncovered stent
- Radiopaque markers for Endoscopic guidance
- Side by Side stenting (2 introducers in one working channel of Endoscope)
- 0.025" Guidewire / 5.9Fr. Delivery device

Endoscopic Application

| | Stent(mm) | | | Delivery device | |
|---|-----------|----------------|---------------|-----------------|--------------|
| | Diameter | Usable length* | Total length* | Length(mm) | Diameter(Fr) |
| BNHS-06-040-180 ⋮ BNHS-06-120-180 | 6 | - | 40~120* | 1800 | 5.9 |
| BNHS-08-040-180 ⋮ BNHS-08-120-180 | 8 | - | 40~120* | 1800 | 5.9 |

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device

Distal Proximal

(NNN) Non covered
(CCC) Fully covered
(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length
② Stent Diameter ④ Delivery Device Length

HANAROSTENT® Benefit™ Biliary (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Uncovered stent
- Radiopaque markers for Percutaneous guidance
- Side by Side stenting (2 introducers in one working channel of Endoscope)
- 0.025" Guidewire / 5.9Fr. Delivery device

Percutaneous Application

| | Stent(mm) | | | Delivery device | |
|---|-----------|----------------|---------------|-----------------|--------------|
| | Diameter | Usable length* | Total length* | Length(mm) | Diameter(Fr) |
| BNHS-06-040-060 ⋮ BNHS-06-120-060 | 6 | - | 40~120* | 600 | 5.9 |
| BNHS-08-040-060 ⋮ BNHS-08-120-060 | 8 | - | 40~120* | 600 | 5.9 |

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device

Distal Proximal

(NNN) Non covered
(CCC) Fully covered
(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length
② Stent Diameter ④ Delivery Device Length