


HANAROSTENT® Esophagus Upper (CCC)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :


- Flare head on the proximal to be positioned near to Upper Esophageal Sphincter
- Fully covered and lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
ECP18060-X070 ⋮ ECP18170-X070	24-18-26	32~142*	60~170*	700	6/18
ECP20060-X070 ⋮ ECP20170-X070	26-20-28	32~142*	60~170*	700	6/18

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device



Proximal

(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length