


HANAROSTENT® Esophagus Asymmetric (CCC)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :


- Shorter dumbbell head on the proximal to be positioned near to Upper Esophageal Sphincter
- Fully covered and lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
ECN18060-X070 ⋮ ECN18170-X070	24-18-24	30~140*	60~170*	700	6/18
ECN20060-X070 ⋮ ECN20170-X070	26-20-26	30~140*	60~170*	700	6/18
ECN24060-Z070 ⋮ ECN24170-Z070	30-24-30	30~140*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④
① Product Name ③ Stent Length
② Stent Diameter ④ Delivery Device Length