


HANAROSTENT® Biliary (NC)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Large flange on the proximal part for Anti-Migration toward the liver
- Uncovered on the distal part for Anti-Migration toward the stomach
- Radiopaque markers for fluoroscopic guidance
- Lasso on proximal end

Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BPD-10-060-180 ⋮ BPD-10-100-180	10-20	54~94*	60~100*	1800	8.5

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~100 : 60,70,80,90,100)

• Based on Delivery Device

Distal Proximal

(NNN) Non covered
(CCC) Fully covered
(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length
② Stent Diameter ④ Delivery Device Length