


HANAROSTENT®

Plumber™



Indication :

- Transgastric or transduodenal endoscopic drainage of symptomatic pancreatic pseudocysts and symptomatic Walled Off Necrosis.

Features :

- Large wheels for Anti-Migration
- Fully covered stent for prevent leakage
- Optimal radial force of stent body
- Radiopaque markers for fluoroscopic guidance

Endoscopic Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BCF-10-020-180 ⋮ BCF-10-040-180	22-10-22	13~33*	20~40*	1800	10.5
BCF-12-020-180 ⋮ BCF-12-040-180	24-12-24	13~33*	20~40*	1800	10.5
BCF-14-020-180 ⋮ BCF-14-040-180	26-14-26	12~32*	20~40*	1800	10.5
BCF-16-020-180 ⋮ BCF-16-040-180	28-16-28	12~32*	20~40*	1800	10.5

※ Sizing and availability varies by country

* Increment of 10mm (Total length 20~40 : 20,30,40)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length