


CHOOSTENT® Esophagus (CCC)



Indication :

- Palliative treatment of esophageal stricture and/or trachea-esophageal fistula caused by malignant tumors

Features :

- Dumbbell heads on the both ends for Anti-Migration
- To be positioned in the Lower Esophageal Sphincter
- Fully covered and lassos on the both ends
- Segmented body structure to conform to peristaltic movement

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
NES-18-060-070 ⋮ NES-18-170-070	24-18-24	20~130*	60~170*	700	6/18
NES-20-060-070 ⋮ NES-20-170-070	26-20-26	20~130*	60~170*	700	6/18
NES-22-060-070 ⋮ NES-22-170-070	28-22-28	20~130*	60~170*	700	8/24
NES-24-060-070 ⋮ NES-24-170-070	30-24-30	20~130*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length