


HANAROSTENT® Colon/Rectum (CCC)



Indication :

- Palliative treatment of colo-rectal strictures caused by malignant tumors

Features :

- Dumbbell heads on the both ends for Anti-Migration
- To be positioned in the descending colon
- Fully covered and lassos on the both ends

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
CCB22060-Z070 ⋮ CCB22170-Z070	30-22-30	20~130*	60~170*	700	8/24
CCB24060-Z070 ⋮ CCB24170-Z070	32-24-32	20~130*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device

Distal Proximal

(NNN) Non covered
(CCC) Fully covered
(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name ③ Stent Length
② Stent Diameter ④ Delivery Device Length