


HANAROSTENT® Colon/Rectum TTS (CCC)



Indication :

- Palliative treatment of colo-rectal strictures caused by malignant tumors

Features :

- Flare heads on the both ends for Anti-Migration
- To be positioned in the colon and rectum
- Fully covered and lassos on the proximal end
- Through the Scope (TTS)

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
CCBA-20-060-230 ⋮ CCBA-20-140-230	26-20-26	24~104*	60~140*	2300	3.4/10.2

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~140 : 60,70,80,90,100,110,120,130,140)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length