


CHOOSTENT® Colon/Rectum Asymmetric (CCC)



Indication :

- Palliative treatment of colo-rectal strictures caused by malignant tumors

Features :

- Shorter head on the proximal to be position close to rectum
- Fully covered and lassos on the both ends
- Segmented body structure to conform to peristaltic movement

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(mm/Fr)
CCI22060-Z070 ⋮ CCI22170-Z070	30-22-30	30~140*	60~170*	700	8/24
CCI24060-Z070 ⋮ CCI24170-Z070	32-24-32	30~140*	60~170*	700	8/24

※ Sizing and availability varies by country

* Increment of 10mm (Total length 60~170 : 60,70,80,90,100,110,120,130,140,150,160,170)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length