


HANAROSTENT® Biliary (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Flare on the both ends for Anti-Migration
- Uncovered stent
- Radiopaque markers for fluoroscopic guidance
- Flexible body structure and High Radial Force

Endoscopic Application


	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BNA08040-G180 ⋮ BNA08120-G180	11.5-8-11.5	26~106*	40~120*	1800	8
BNA10040-G180 ⋮ BNA10120-G180	13.5-10-13.5	26~106*	40~120*	1800	8

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device		(NNN) Non covered	XXX-00-000-000	XXX 00 000-X000
		(CCC) Fully covered	① ② ③ ④	① ② ③ ④
		(NCN) Partially covered	① Product Name	③ Stent Length
			② Stent Diameter	④ Delivery Device Length

HANAROSTENT® Biliary (NNN)



Indication :

- Palliative treatment of bile duct stricture caused by malignant tumors

Features :

- Flare on the both ends for Anti-Migration
- Uncovered stent
- Radiopaque markers for fluoroscopic guidance
- Flexible body structure and High Radial Force

Percutaneous Application

	Stent(mm)			Delivery device	
	Diameter	Usable length*	Total length*	Length(mm)	Diameter(Fr)
BNA08040-H060 ⋮ BNA08120-H060	11.5-8-11.5	26~106*	40~120*	600	8
BNA10040-H060 ⋮ BNA10120-H060	13.5-10-13.5	26~106*	40~120*	600	8

※ Sizing and availability varies by country

* Increment of 10mm (Total length 40~120 : 40,50,60,70,80,90,100,110,120)

• Based on Delivery Device



(NNN) Non covered

(CCC) Fully covered

(NCN) Partially covered

XXX-00-000-000 XXX 00 000-X000

① ② ③ ④ ① ② ③ ④

① Product Name

② Stent Diameter

③ Stent Length

④ Delivery Device Length